EMR-validated chair scale with precise graduation



seca 954

- + 0.1 pound graduation for exact measurement while sitting
- + Fold-up arm and footrests for easy access
- + Ergonomic seat for increased patient comfort
- + Suitable for patients affected by obesity with up to 660 pounds weighing capacity
- + Sophisticated engineering and robust construction for a long lifespan
- + Parking brake for secure position
- + EMR-validated: transmit measured data directly to an electronic medical record system or printer



EMR-validated chair scale with precise graduation



Highly accurate and convenient

With its fine graduation, the seca 954 is undisputed king of the chair scales. However, it is not only the most accurate, but also belongs to the premium class in terms of comfort. Its low weight, maneuverability and folding foot and armrests make daily work considerably easier. The chair scale is designed so that it is particularly quick and easy to clean. This means that all parts of the seca 954 can be reached by hand, and there are no edges, grooves or other traps for pathogens.

Capacity	660 lbs / 300 kg
Graduation	0.1 lb < 330 lbs > 0.2 lb / 50 g < 150 kg > 100 g
Dimensions (WxHxD)	22.2 x 35.3 x 38.5 in / 563 x 897 x 978 mm
Net weight	54 lbs / 24.5 kg
Power supply	Power adapter, Rechargeable batteries
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, RESET, Automatic switch-off, lbs/kg switch-over, Adjustable damping, BMI
Interfaces	Wi-Fi / Ethernet optional with seca 452
Accessories	+ Interface module seca 452 + Rechargeable battery seca 454



Why you can rely on seca quality

The seca 954 chair scale has an extremely stable round tubular frame that has been processed to a high level of craftsmanship. Due to the powder coating, the hardware is protected against knocks and scratches. The rear swivel rollers are firmly anchored in the tube and countered. In addition, the rollers are conductive, which prevents electrostatic charging. As before, the chair scales and all other multifunctional scales are manufactured, tested and sent across the world at seca headquarters in Hamburg. In this way, we guarantee the functionality, durability and robustness of our products first-hand.



seca integration solutions: How to significantly improve workflows

EMR-validation means that the seca 954 can be integrated into the hospital information system via the network. In this way, measured data can be sent directly to the digital patient file (EMR) at the touch of a button, where they are immediately available to all users. With an additionally connected USB barcode scanner, the patient and user are identified and the measurement results are always correctly assigned. This saves you time, increases process reliability and eliminates transmission errors.

More information on seca integration solutions can be found at: integration.seca.com